



# Indigo® System

Mechanical Thrombectomy

## CAT RX

As part of the Indigo Aspiration System, the Indigo CAT RX Aspiration Catheters and Indigo Separator™ 4 are indicated for the removal of fresh, soft emboli and thrombi from vessels in the coronary and peripheral vasculature.

Penumbra ENGINE™



# CAT RX

GET THE CLOT OUT FAST<sup>a</sup>

Advanced Trackability

Seven Material Transitions

Enhanced Deliverability

Proximal Laser Cut Hypotube

Maximized Clot Removal

Large Lumen | 6 F Guide Compatibility

Rapid Exchange


0.014" Guidewire Compatibility

## DATA PRESENTED AT ACC 2019

Study	“Initial Experience with a Mechanical Aspiration Catheter for Thrombus Removal During Percutaneous Intervention: A Multicenter Retrospective Case Series”		
Key Results	<div><div>93.2%</div><div>Post-CAT RX TIMI 3 flow</div></div>	<div><div><div>35 SEC</div></div><div>Median aspiration time with CAT RX</div></div>	<div><div>0%</div><div>Incidence of stroke</div></div>
Design	<div><div>59 patients</div><div>With fresh, soft thrombi in the RCA, LAD, or Circumflex</div></div>	<div><div>76.3%</div><div>Of patients with an occlusion of TIMI 0 flow pre-procedure</div></div>	
Conclusion	Results are encouraging and can be used to plan future prospective trials.		


Mathews SJ, Brown C, Kolski B, et al. Initial experience with a mechanical aspiration catheter for thrombus removal during percutaneous intervention: a multicenter retrospective case series. Poster presented at: American College of Cardiology’s 68th Annual Scientific Session; March 16-18, 2019; New Orleans, LA.

### Removal of Thrombus from Right Coronary Artery





Dr. Suhail Dohad, Cedars-Sinai Medical Center, CA

### Removal of Thrombus from Left Anterior Descending Artery



Dr. Jay Mathews, Manatee Memorial Hospital, FL

## Mechanical Power Aspiration with Penumbra ENGINE<sup>b</sup>




<sup>b</sup>. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance.

<sup>a</sup>. Time represents duration of angiogram. Duration of angiogram provided by consenting physician. Images used with permission. Consent on file at Penumbra, Inc. Individual results may vary depending on a variety of patient-specific attributes.




## Removal of Thrombus from Right Coronary Artery




Dr. Nainesh Patel, Lehigh Valley Health Network, PA

## Removal of Thrombus from LIMA Graft




Dr. Larry Diaz, Metro Health-University of Michigan Health, MI

## Removal of Thrombus from Left Anterior Descending Artery



Dr. Branavan Umakanthan, Nevada Heart and Vascular Center, NV

## Removal of Thrombus from Circumflex Artery



Dr. Brian Kolski, St. Joseph Hospital, CA


## Here's what physicians are saying about CAT RX:




Jon George  
@jgeorgemd



To aspirate or not? Still doing it in my acute MI patients with large thrombotic debris. You be the judge. [@GPAngioClub](#) [@ACVCPhila](#) [@penumbrainc](#)




8 July 2019



Ethan Korngold, MD  
@ekgpdx



CatRX via 6F guide. Surprisingly deliverable and the continuous suction is easier to use than manual aspiration devices. I will use if extraordinary thrombus burden. [@bkolskk](#)




19 Jun 2019



Shariq Shamim, MD, FACC  
@ShariqShamimMD



you have to use the device to believe it. You can see thrombus in connection line or canister rest is filled with blood. Really good aspiration when large clot. Cc: [@bkolskk](#) [@AntoniousAttall](#)



18 Aug 2019

Connect with us on Twitter

[@PenVascular](#)

a. Time represents duration of angiogram. Duration of angiogram provided by consenting physician. Images used with permission. Consents on file at Penumbra, Inc. Individual results may vary depending on a variety of patient-specific attributes.

# Ordering Information

Indigo® Catheter Kits							
Catalog Number	Description	Proximal OD (F)	Distal OD (F)	Compatibility (F) (Sheath or Guide)	Working Length (cm)	Wire Platform (in.)	Compatible Penumbra Devices
CATRXXKIT	Indigo CAT RX + Large Lumen Aspiration Tubing	–	5.3	6.0 Sheath or Guide	140	.014	Separator™ 4
CAT8XTORQ115KIT	Indigo 8 XTORQ Tip + Dynamic Aspiration Tubing	8.0	8.0	8.0 F Sheath	115	.014–.038	Separator 8
CAT8TORQ85KIT	Indigo 8 TORQ Tip + Dynamic Aspiration Tubing	8.0	8.0	8.0 F Sheath	85	.014–.038	Separator 8
CAT8STR85KIT	Indigo 8 Straight Tip + Dynamic Aspiration Tubing	8.0	8.0	8.0 F Sheath	85	.014–.038	Separator 8
CAT6KIT	Indigo 6 + Dynamic Aspiration Tubing	6.0	6.0	6.0 F Sheath	135	.014–.038	Separator 6
CAT5KIT	Indigo 5 + Dynamic Aspiration Tubing	6.0	5.0	6.0 F Sheath	132	.014–.038	Separator 5
CAT3KIT	Indigo 3 + Dynamic Aspiration Tubing	4.1	3.4	5.0 F Sheath	150	.014–.025	Separator 3
CATD	Indigo D + Large Lumen Aspiration Tubing	8.0	8.0	8.0 F Sheath	50	.014–.038	Separator D

Indigo Separators					Accessories		
Catalog Number	Description	Distal OD (in.)	Total Length (cm)	Compatible Penumbra Devices	Catalog Number	Description	Compatible Penumbra Devices
SEPC4	Separator 4	.035	200	CAT RX	PMXENGN	Penumbra ENGINE™	Penumbra ENGINE Canister
SEP8	Separator 8	.072	150	CAT8	IAPS3	Penumbra ENGINE Canister	Penumbra ENGINE
SEP6	Separator 6	.055	175	CAT6			
SEP5	Separator 5	.045	175	CAT5			
SEP3	Separator 3	.028	190	CAT3			
SEPD	Separator D	.072	90	CATD			

## INDIGO Aspiration Catheters and SEPARATORS – Indication For Use

As part of the INDIGO Aspiration System, the INDIGO Aspiration Catheters and SEPARATORS are indicated for the removal of fresh, soft emboli and thrombi from vessels of the peripheral arterial and venous systems.

### INDIGO Aspiration Tubing – Indication For Use

As part of the INDIGO Aspiration System, the INDIGO Sterile Aspiration Tubing is indicated to connect the INDIGO Aspiration Catheters to the Penumbra Aspiration Pump.

### Penumbra Aspiration Pump – Indication For Use

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

### Contraindications

Not for use in the coronaries or the neurovasculature.

### Warnings

- The INDIGO Aspiration System should only be used by physicians who have received appropriate training in interventional techniques.
- Do not advance, retract or use any component of the INDIGO System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter or SEPARATOR against resistance may result in damage to the device or vessel.
- Do not use the INDIGO Aspiration System with a pump other than the Penumbra Aspiration Pump.

### Precautions

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the "Use By" date.
- Use the INDIGO Aspiration System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the INDIGO Aspiration Tubing valve is open for only the minimum time needed to remove thrombus. Excessive aspiration or failure to close the INDIGO Aspiration Tubing valve when aspiration is complete is not recommended.

- The INDIGO SEPARATOR is not intended for use as a guidewire. If repositioning of the INDIGO Aspiration Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate guidewire using standard microcatheter and guidewire techniques.
- Do not use automated high-pressure contrast injection equipment with the INDIGO Aspiration Catheter because it may damage the device.

### Potential Adverse Events

Possible complications include, but are not limited to, the following: allergic reaction and anaphylaxis from contrast media; acute occlusion; air embolism; arteriovenous fistula; death; device malfunction; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; inability to completely remove thrombus; infection; hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; vessel spasm, thrombosis, dissection, or perforation; intimal disruption; myocardial infarction; emergent surgery; fibrillation; hypotension; respiratory failure; peripheral thromboembolic events.

### INDIGO CAT RX Aspiration Catheters and

#### INDIGO SEPARATOR 4 – Indication For Use

As part of the INDIGO Aspiration System, the INDIGO CAT RX Aspiration Catheters and INDIGO SEPARATOR 4 are indicated for the removal of fresh, soft emboli and thrombi from vessels in the coronary and peripheral vasculature.

### INDIGO Aspiration Tubing – Indication For Use

As part of the INDIGO Aspiration System, the INDIGO Sterile Aspiration Tubing is indicated to connect the INDIGO CAT RX Aspiration Catheters to the Penumbra Aspiration Pump.

### Penumbra Aspiration Pump – Indication For Use

The Penumbra Aspiration Pump is indicated as a vacuum source for Penumbra Aspiration Systems.

### Contraindications

- The INDIGO Aspiration System is contraindicated in:
  - The removal of fibrous, adherent or calcified material (e.g. chronic clot, atherosclerotic plaque)
  - The cerebral vasculature

### Warnings

- The INDIGO Aspiration System should only be used by physicians who have received appropriate training in interventional techniques.
- Do not advance, retract, or use any component of the INDIGO Aspiration System against resistance without careful

assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torquing or forced insertion of the catheter or SEPARATOR against resistance may result in damage to the device or vessel.

- Do not use the INDIGO Aspiration System with a pump other than the Penumbra Aspiration Pump.

### Precautions

- The safety and effectiveness of this device for use in the treatment of ST-Elevation Myocardial Infarction (STEMI) has not been established. Complications from the use of this device in this manner could lead to death, permanent impairment, and/or the need for emergency medical intervention.**

- The device is intended for single use only. Do not resterilize or reuse. Resterilization and/or reuse may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target vasculature location.
- Do not use kinked or damaged devices. Do not use open or damaged packages. Return all damaged devices and packaging to the manufacturer/distributor.
- Use prior to the "Use By" date.
- Use the INDIGO Aspiration System in conjunction with fluoroscopic visualization.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the INDIGO Aspiration Tubing valve is open for only the minimum time needed to remove thrombus. Excessive aspiration or failure to close the INDIGO Aspiration Tubing valve when aspiration is complete is not recommended.
- The INDIGO SEPARATOR 4 is not intended for use as a guidewire.
- If repositioning of the INDIGO CAT RX Aspiration Catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate guidewire using standard guidewire techniques.
- Do not use INDIGO SEPARATOR 4 to macerate or retrieve thrombus distal to the catheter tip. INDIGO SEPARATOR 4 is intended to be used with INDIGO CAT RX Aspiration Catheter to clear the distal end of the catheter lumen should it be blocked with thrombus.
- Do not use automated high-pressure contrast injection equipment with the INDIGO CAT RX Aspiration Catheter because it may damage the device.

### Potential Adverse Events

Possible complications include, but are not limited to, the following: allergic reaction and anaphylaxis from contrast media; acute occlusion; air embolism; arteriovenous fistula; death; device malfunction; distal embolization; emboli; false aneurysm formation; hematoma or hemorrhage at access site; inability to completely remove thrombus; infection; hemorrhage; ischemia; kidney damage from contrast media; neurological deficits including stroke; vessel spasm,

thrombosis, dissection, or perforation; intimal disruption; myocardial infarction; emergent surgery; fibrillation; hypotension; respiratory failure; peripheral thromboembolic events.

### PENUMBRA ENGINE – Indication For Use

The PENUMBRA ENGINE is indicated as a vacuum source for Penumbra Aspiration Systems.

### Contraindications

There are no contraindications.

### Warnings/Precautions

- The canister is intended for single use only. Do not reuse. Re-use may result in canister cracking or vacuum filter blockages, which may result in the inability to aspirate.
- Do not block bottom air vents. Unit may overheat and shut off or fail to restart if run for extended periods of time without airflow.
- To avoid the risk of electrical shock, this equipment must only be connected to a supply mains with protective earth.
- Do not position the PENUMBRA ENGINE so that it is difficult to remove the power cord. The means of mains disconnect is to remove the power cord.
- Only use replacement fuse with correct rating (see Table 1 for fuse rating).
- Remove and service the PENUMBRA ENGINE if liquids or solids have been drawn into the PENUMBRA ENGINE.
- Do not use in the presence of a flammable anesthetic mixture with air or nitrous oxide.
- Do not use in an oxygen rich environment.
- To prevent fire or shock hazard, use a replacement power cord of equal rating.
- Do not re-infuse blood or fluid from the canister back into the patient.
- Do not use petroleum based compounds, acids, caustics, or chlorinated solvents to clean or lubricate any parts. It will reduce the service life of the PENUMBRA ENGINE. Use only water-based solvents for cleaning.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the PENUMBRA ENGINE. Otherwise, this could result in degradation of the performance of this equipment.
- Common emitters (such as RFID emitters, security systems, diathermy equipment, and portable transmitters) should not be used in close proximity to the PENUMBRA ENGINE as they can interfere with and result in degradation of the performance of the equipment.
- Equipment is not safe for MR use.
- No modification of this equipment is allowed.

**Penumbra** 

www.penumbrainc.com

**Penumbra, Inc. USA**  
One Penumbra Place  
Alameda, CA 94502  
USA  
T 1.888.272.4606  
T 1.510.748.3200  
F 1.510.748.3232  
order@penumbrainc.com  
info@penumbrainc.com

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance. Images used with permission. Consents on file at Penumbra, Inc. Individual results may vary depending on a variety of patient-specific attributes. Please contact your local Penumbra representative for more information.

Copyright ©2019 Penumbra, Inc. All rights reserved. The Penumbra P logos, Indigo, Separator, and Penumbra ENGINE are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. All other trademarks are the property of their respective owners. 15996, Rev. A 08/19 USA